OnGuard[®] 2 Closed System Transfer Device (CSTD)
Containing Chemfort[™] CSTD Reconstitution and Administration Components

P-8749-2 REV. 09/22

VIAL ADAPTORS (VA) AND CONVERTOR 13 mm

DRUG

DILUENT/SALINE

13 mm	20 mm	28 mm
VA 20 mm + convertor 13 mm	VA 20 mm	VA 28 mm

SYRINGE ADAPTOR (SA)

SYRINGE ADAPTOR LOCK (SAL)

LUER LOCK ADAPTOR (LLA)

BAG ADAPTOR SP (BASP)

A

VIAL ADAPTOR CONNECTION
Each Vial Adaptor fits a specific vial neck

2a

20mm VA / 28mm VA

3a

CLICK!

2b

13mm

3b

CLICK!

4b

CLICK!

20mm / 28mm

Connect VA spike
Once connected, do not disconnect

13mm

Invert Vial; Push into 13mm convertor;
Connect 20mm VA spike
Once connected, do not disconnect

B

SA/SAL CONNECTION

2a

Connect SA to syringe
Once connected, do not disconnect

2b

Connect SAL to syringe

C

BAG ADAPTOR SP

2

Insert BASP to Infusion Bag
Once connected, do not disconnect

3

D

DRUG RECONSTITUTION (powder for injection)

1a

13mm / 20mm / 28mm
Remove caps

2a

Connect SA/SAL to VA of diluent, invert and withdraw diluent

1b

Connect SA/SAL to BASP; Withdraw diluent; Pinch wings to disconnect

2b

3b

4

Connect SA/SAL of diluent to VA of lyophilized vial, inject diluent; dissolve drug

5

6

7

E

DRUG TRANSFER

1

13mm / 20mm / 28mm
Remove caps

2

Connect SA/SAL to the VA of drug and invert to withdraw

3

Withdraw drug while vial is upwards
Discard air bubbles. Do not inject fluid upwards

4

Pinch wings on SA/SAL to disconnect

5

Re-Cap SA/SAL and VA for future use

F

DRUG INFUSION TO IV BAG

1

Connect SA/SAL of drug to BASP and inject drug

2

3

Flush with 2-3 mL of solution

4

5

Pinch wings on SA/SAL to disconnect

G

ADMINISTRATION-IV DRIP

1a

Open Spike Port tail and insert primed Set Spike
Once connected, do not disconnect

2a

H

ADMINISTRATION- IV PUSH/BOLUS

1

Connect LLA to IV Female Luer Lock (FLL) port
Once connected, do not disconnect

2

3

Connect SA/SAL of Drug to LLA and inject

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Explanation of symbols taken from ISO 15223-1:2021

ref: 5.1.1		Manufacturer	ref: 5.1.6		Catalog number	ref: 5.2.8		Do not use if package damaged	ref: 5.3.2		Keep away from sunlight	ref: 5.4.3		Consult instructions for use
ref: 5.1.4		Use-by date	ref: 5.2.3		Sterilized using ethylene oxide	ref: 5.2.11		Single sterile barrier system	ref: 5.3.4		Keep dry	ref: 5.6.3		Non-pyrogenic
ref: 5.1.5		Batch code	ref: 5.2.6		Do not resterilize	ref: 5.2.14		Single sterile barrier system with protective packaging outside	ref: 5.4.2		Do not re-use	ref: 5.7.8	 21 Yitzhak Sadeh St. Raanana Israel 4352621 Global Translations	Translation